


Intensiv diamond metal strips with handles

Manual contouring, finishing and polishing of proximal surfaces




- Contouring, finishing and polishing in just one process.
- Short strip for improved finger support
- Sterilizable and reusable



Finishing and polishing
Interproximal anterior restorations



REF PXA4015/6



Finishing and polishing
All interproximal anterior restorations




REF PX4015/6



Morphological corrections



REF PXCC8040/6



Gloss polishing



REF PXPO1508/6


Intensiv Metal Diamond Strips

By phone : **+ 41 91 986 50 50** or by e-mail: **customer.care@intensiv.ch**

Product	REF	Price	Quantity
			
Intensiv Diamond Strips Assortment, Trial Kit Content: Intensiv ProxoStrip Anterior, 1 piece Intensiv ProxoStrip, 2 pieces Intensiv ProxoContour Coarse, 2 pieces Intensiv ProxoPolish, 1 piece	PXT/6		
			
Intensiv ProxoStrip Anterior	PXA4015/6		
Intensiv ProxoStrip Anterior	PXA4015/12		
			
Intensiv ProxoStrip	PX4015/6		
Intensiv ProxoStrip	PX4015/12		
			
Intensiv ProxoContour Coarse	PXCC8040/6		
			
Intensiv ProxoPolish	PXPO1508/6		

My dental clinic

Date _____ Signature _____



☐ Please send an Intensiv Catalogue

AUTOMATIC E-MAIL ENQUIRY

Intensiv SA
 Via al Molino 107
 6926 Montagnola
 Switzerland

Tel. + 41 91 986 50 50
 order@intensiv.ch
 www.intensiv.ch

